

MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit

High throughput isolation of viral nucleic acid (RNA and DNA) from biofluids and transport media

Catalog Number A42352

Pub. No. MAN0018073 Rev. B.0

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from thermofisher.com/support.

Product description

The Applied Biosystems™ MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit is developed for scalable, rapid purification of high-quality nucleic acid (RNA and DNA) from virus and easy to lyse bacteria in biofluid and transport media samples. You can use the nucleic acid purified with this kit in a broad range of molecular biology downstream applications, such as sequencing and qPCR. This protocol guides users through automated isolation of nucleic acid using the KingFisher™ Flex and the KingFisher™ Duo Prime instruments.

Contents and storage

Reagents that are provided in the kit are sufficient for 100 reactions with standard volume input or 20 reactions with large volume input.

Table 1 Components of MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit (Cat. No. A42352)

Component	Amount	Storage
Binding Solution	53 mL	15°C to 25°C
Wash Buffer	100 mL	
Elution Solution	10 mL	
Proteinase K	1 mL	
Total Nucleic Acid Binding Beads	2 mL	

For 1,000 reaction volume, use Cat. No. A42359 (Binding Solution), A42360 (Wash Buffer), A42364 (Elution Solution), A42363 (Proteinase K), A42362 (Binding Beads).

Required materials not supplied

Unless otherwise indicated, all materials are available through thermofisher.com. MLS: Fisher Scientific (fisherscientific.com) or other major laboratory supplier.

Item	Source
Instrument	
Magnetic particle processor (one of the following, depending on quantity/volume of sample to be processed):	
<i>For standard volume sample</i> ^[1] : KingFisher™ Flex Magnetic Particle Processor with 96 Deep-Well Head	5400630
<i>For large volume sample</i> ^[2] : KingFisher™ Flex Magnetic Particle Processor with 24 Deep-Well Head	5400640
KingFisher™ Duo Prime Magnetic Particle Processor	5400110
Consumables	
Deep-well plates:	
<i>For standard volume sample</i> ^[1] : KingFisher™ deep-well 96 plate	95040450
<i>For large volume sample</i> ^[2] : KingFisher™ Flex 24 deep-well plate	95040470
96-well standard plates (for use with KingFisher™ Flex only; tip comb placement and eluate storage):	
KingFisher™ 96 KF plate	97002540
Tip comb, compatible with the magnetic particle processor used:	
KingFisher™ Duo Prime 12-tip comb, for use with KingFisher™ deep-well 96 plate	97003500
KingFisher™ Duo Prime 6-tip comb, for use with KingFisher™ Flex 24 deep-well plate	97003510
KingFisher™ 96 tip comb for deep-well magnets, KingFisher™ Flex protocol only	97002534
KingFisher™ Flex 24 deep-well tip comb and plate, KingFisher™ Flex protocol only	97002610
Elution strip (for use with KingFisher™ Duo Prime only; elution step):	
KingFisher Duo elution strip	97003520
KingFisher Duo cap for elution strip	97003540

Item	Source
Equipment	
Adjustable micropipettors	MLS
Multi-channel micropipettors	MLS
Materials	
MicroAmp™ Clear Adhesive Film	4306311
Conical Tubes (15 mL)	AM12500
Conical Tubes (50 mL)	AM12501
Reagent reservoirs	MLS
Nonstick, RNase-Free Microfuge Tubes, 1.5 mL	AM12450
Nonstick, RNase-Free Microfuge Tubes, 2.0 mL	AM12475
Reagents	
Ethanol, 100% (molecular biology grade)	MLS
Nuclease-free Water	AM9932

^[1] Standard volume sample is 200–400 µL.

^[2] Large volume sample is 500 µL–2 mL.

General guidelines

- Perform all steps at room temperature (20–25°C), unless otherwise noted.
- Precipitates can occur if the Binding Solution is stored when room temperature is too cold. If there are precipitates, warm the Binding Solution at 37°C and gently mix to dissolve the precipitates. Avoid creating bubbles.

- Reagent Mix tables are sufficient for a single reaction. To calculate volumes for other sample numbers, see the per-well volume and add at least a 10% overage.
- (Optional): To prevent evaporation and contamination, cover the prepared processing plates with paraffin film or MicroAmp™ Clear Adhesive Film until they are loaded into the instrument.

Guidelines for Binding Bead Mix

- Vortex Binding Beads thoroughly before each use.
- Ensure that the beads stay fully mixed within the solution during pipetting.
- Avoid creating bubbles during mixing and aliquoting.
- Binding/Bead Mix is very viscous so pipet with care to ensure that the correct volume is added to the sample.

Before first use of the kit

- Prepare 80% Ethanol from 100% absolute Ethanol and Nuclease-Free Water.
 - For standard volume input: Prepare enough for 1.5 mL per reaction.
 - For large volume input: Prepare enough for 6 mL per reaction.

Perform total nucleic acid purification using KingFisher™ Flex (standard volume: 200–400 µL)

1 Set up the instrument

- a. Ensure that the instrument is set up with the proper magnetic head and the proper heat block, as indicated in the following table.

Component	Type
Magnetic head	96 deep-well magnetic head
Heat block	96 deep-well heat block

IMPORTANT! Failure to use the proper magnetic head and heat block results in lower yields and potential harm to the instrument.

- b. Ensure that the proper program (**MVP_Flex**) has been downloaded from the product page and loaded onto the instrument.

2 Set up the processing plates

Set up the Wash, Elution, and Tip Comb Plates outside the instrument according to the following table.

Plate ID	Plate position	Plate type	Reagent	Volume per well
Wash 1 Plate	2	Deep-well	Wash Buffer	1,000 µL
Wash 2 Plate	3	Deep-well	80% Ethanol	1,000 µL
Wash 3 Plate	4	Deep-well	80% Ethanol	500 µL
Elution Plate	5	Deep-well	Elution Solution	50–100 µL
Tip Comb	6	Place a 96 Deep-well Tip Comb in a Standard Plate		

3

Prepare Binding Bead Mix

- a. Vortex Beads vigorously to ensure they are homogenous.
- b. Prepare Binding Bead Mix according to the following table and sample input volume:

Component	Volume per well ^[1]
Binding Solution	530 µL
Total Nucleic Acid Magnetic Beads	20 µL
Total volume	550 µL

^[1] Use 10% Overage calculation when making a master mix for use with multiple samples.

- c. Mix well by inversion, then store at room temperature.

4

Digest with Proteinase K, then elute nucleic acid

- a. Add 10 µL of Proteinase K to each well of a Deep-well 96-well plate.
This plate is the Sample Plate.
- b. Add 200–400 µL of each sample to wells with Proteinase K in the Sample Plate.
Note: Recommend up to 200 µL input for whole blood.
- c. Invert Binding Bead Mix gently to mix, then add 550 µL to each sample in the Sample Plate.
Note: Remix the Binding Bead Mix by inversion frequently during pipetting to ensure even distribution of beads to all samples or wells. The mixture containing the Binding Beads is viscous. Therefore, pipet slowly to ensure that the correct amount is added. DO NOT use a repeat pipet to add to the samples as the high viscosity will cause variations in volume added.
- d. Select the program **MVP_Flex** on the instrument.
- e. Start the run, then load the prepared plates into position when prompted by the instrument.
- f. After the protocol is complete (~25 minutes after start), immediately remove the elution plate from the instrument and cover the plate or transfer the eluate to a tube or plate of choice for final storage.

The purified nucleic acid is ready for immediate use. Alternatively, store the plate at –20°C for long-term storage.

Perform total nucleic acid purification using KingFisher™ Flex (large volume: 500 µL to 2 mL)

1

Set up the instrument

- a. Ensure that the instrument is set up with the proper magnetic head and the proper heat block, as indicated in the following table.

Component	Type
Magnetic head	24 deep-well magnetic head
Heat block	24 deep-well heat block

IMPORTANT! Failure to use the proper magnetic head and heat block results in lower yields and potential harm to the instrument.

- b. Ensure that the proper program (**MVP_Flex_LV**) has been downloaded from the product page and loaded onto the instrument.

2 Set up the processing plates

Set up the Wash, Elution, and Tip Comb Plates outside the instrument according to the following table.

Plate ID	Plate position	Plate type	Reagent	Volume per well
500 µL sample input				
Wash 1 Plate	2	24 Deep-well	Wash Buffer	2,000 µL
Wash 2 Plate	3	24 Deep-well	80% Ethanol	2,000 µL
Wash 3 Plate	4	24 Deep-well	80% Ethanol	2,000 µL
Elution Plate	5	24 Deep-well	Elution Solution	150 µL
Tip Comb	6	Place a 24 Deep-well Tip Comb in a 24 deep-well Plate		
>500 µL – 1 mL sample input				
Wash 1 Plate	2	24 Deep-well	Wash Buffer	4,000 µL
Wash 2 Plate	3	24 Deep-well	80% Ethanol	4,000 µL
Wash 3 Plate	4	24 Deep-well	80% Ethanol	2,000 µL
Elution Plate	5	24 Deep-well	Elution Solution	200 µL
Tip Comb	6	Place a 24 Deep-well Tip Comb in a 24 deep-well Plate		
>1 mL – 2 mL sample input				
Wash 1 Plate	2	24 Deep-well	Wash Buffer	4,000 µL
Wash 2 Plate	3	24 Deep-well	80% Ethanol	4,000 µL
Wash 3 Plate	4	24 Deep-well	80% Ethanol	2,000 µL
Elution Plate	5	24 Deep-well	Elution Solution	250 µL
Tip Comb	6	Place a 24 Deep-well Tip Comb in a 24 deep-well Plate		

3 Prepare Binding Bead Mix

a. Vortex Beads vigorously to ensure they are homogenous.

b. Prepare Binding Bead Mix according to the following table and sample input volume:

Component	Volume per well ^[1]
500 µL sample input	
Binding Solution	1,250 µL
Total Nucleic Acid Magnetic Beads	50 µL
Total volume	1,300 µL
500 µL – 2 mL sample input	
Binding Solution	2,700 µL
Total Nucleic Acid Magnetic Beads	100 µL
Total volume	2,800 µL

^[1] Use 10% Overage calculation when making a master mix for use with multiple samples.

c. Mix well by inversion, then store at room temperature.

**Digest with Proteinase K,
then elute nucleic acid**

- a. Add Proteinase K to each well in the Sample Plate, according to the following table and sample input volume.

For Sample input volume	Add Proteinase K
500 µL	25 µL
>500 µL – 1 mL	50 µL
>1 mL – 2 mL	100 µL

- b. Add 500 µL to 2 mL of each sample to wells with Proteinase K in the Sample Plate.

Note: Recommend up to 1 mL input for whole blood.

- c. Invert Binding Bead Mix gently to mix, then add the appropriate amount of Binding Bead Mix to each sample in the Sample Plate according to the following table and sample input volume.

For Sample input volume	Add Binding Bead Mix
500 µL	1,300 µL
>500 µL – 2 mL	2,800 µL

Note: Remix the Binding Bead Mix by inversion frequently during pipetting to ensure even distribution of beads to all samples or wells. The mixture containing the Binding Beads is viscous. Therefore, pipet slowly to ensure that the correct amount is added. DO NOT use a repeat pipet to add to the samples as the high viscosity will cause variations in volume added.

- d. Select the program **MVP_Flex_LV** on the instrument.
- e. Start the run, then load the prepared plates into position when prompted by the instrument.
- f. After the protocol is complete (~35 minutes after start), immediately remove the elution plate from the instrument and cover the plate or transfer the eluate to a tube or plate of choice for final storage.

The purified nucleic acid is ready for immediate use. Alternatively, store the plate at –20°C for long-term storage.

Perform total nucleic acid purification using KingFisher™ Duo Prime (standard volume: 200–400 µL)

1 Set up the instrument

- a. Ensure that the instrument is set up with the proper magnetic head and the proper heat block, as indicated in the following table.

Component	Type
Magnetic head	12-tip magnetic head
Heat block	12 well heat strip

IMPORTANT! Failure to use the proper magnetic head and heat block results in lower yields and potential harm to the instrument.

- b. Ensure that the proper program (**MVP_Duo**) has been downloaded from the product page and loaded onto the instrument.

2 Set up the Sample Plate and Elution Strip

Set up the Sample Plate and Elution Strip according to the following tables, respectively.

Table 2 Sample plate

Row ID	Plate Row	Reagent	Volume per well
Sample	A	Sample	Varies
—	B	Empty	
Wash 1	C	Wash Buffer	1,000 µL
—	D	Empty	
Wash 2	E	80% Ethanol	1,000 µL
—	F	Empty	
Wash 3	G	80% Ethanol	500 µL
Tip Comb	H	Tip Comb	

Table 3 Elution strip

Row ID	Plate Row	Reagent	Volume per well
Elution	A	Elution Solution	50–100 µL

3 Prepare Binding Bead Mix

- a. Vortex Beads vigorously to ensure they are homogenous.
- b. Prepare Binding Bead Mix according to the following table and sample input volume:

Component	Volume per well ^[1]
Binding Solution	530 µL
Total Nucleic Acid Magnetic Beads	20 µL
Total volume	550 µL

^[1] Use 10% Overage calculation when making a master mix for use with multiple samples.

- c. Mix well by inversion, then store at room temperature.

4 Digest with Proteinase K, then elute nucleic acid

- a. Add 10 µL of Proteinase K to the wells of Row A of the Sample Plate.
- b. Add 200–400 µL of each sample to wells with Proteinase K in Row A of the Sample Plate.
- Note:** Recommend up to 200 µL input for whole blood.

4 Digest with Proteinase K, then elute nucleic acid
(continued)

- c. Invert Binding Bead Mix gently to mix, then add 550 µL to each sample in Row A of the Sample Plate.
Note: Remix the Binding Bead Mix by inversion frequently during pipetting to ensure even distribution of beads to all samples or wells. The mixture containing the Binding Beads is viscous. Therefore, pipet slowly to ensure that the correct amount is added. DO NOT use a repeat pipet to add to the samples as the high viscosity will cause variations in volume added.
- d. Select the program **MVP_DUO** on the instrument.
- e. Start the run, then load the Elution Strip and Sample Plate into position when prompted by the instrument.
- f. After the protocol is complete (~25 minutes after start), immediately remove the Elution Strip from the instrument. Cover with the Elution Strip Cap for temporary storage, or transfer the eluate to a tube or plate of choice for final storage.

The purified nucleic acid is ready for immediate use. Alternatively, store the plate at -20°C for long-term storage.

Perform total nucleic acid purification using KingFisher™ Duo Prime (large volume: 500 µL – 2 mL)

1 Set up the instrument

- a. Ensure that the instrument is set up with the proper magnetic head and the proper heat block, as indicated in the following table.

Component	Type
Magnetic head	6-tip magnetic head
Heat block	Both 6 well heat strips

IMPORTANT! Failure to use the proper magnetic head and heat block results in lower yields and potential harm to the instrument.

- b. Ensure that the proper program (**MVP_Duo_LV**) has been downloaded from the product page and loaded onto the instrument.

2 Set up the processing plates

Set up the Sample Plate and Elution Plate according to the following tables and sample input volume.

Table 4 24 deep-well plate layout (500 µL sample input)

Row ID	Plate Row	Reagent	Volume per well
Sample plate			
Sample	A	Sample	Varies
Wash 1	B	Wash Buffer	2,000 µL
Wash 2	C	80% Ethanol	2,000 µL
Wash 3	D	80% Ethanol	2,000 µL
Elution plate			
Elution	A	Elution Solution	150 µL
Tip Comb	B	6-Tip Comb in 24 Deep-well plate	

2 Set up the processing plates *(continued)*

Table 5 24 deep-well plate layout (>500 µL – 1 mL sample input)

Row ID	Plate Row	Reagent	Volume per well
Sample plate			
Sample	A	Sample	Varies
Wash 1	B	Wash Buffer	4,000 µL
Wash 2	C	80% Ethanol	4,000 µL
Wash 3	D	80% Ethanol	2,000 µL
Elution plate			
Elution	A	Elution Solution	200 µL
Tip Comb	B	6-Tip Comb in 24 Deep-well plate	

Table 6 24 deep-well plate layout (>1 mL – 2 mL sample input)

Row ID	Plate Row	Reagent	Volume per well
Sample plate			
Sample	A	Sample	Varies
Wash 1	B	Wash Buffer	4,000 µL
Wash 2	C	80% Ethanol	4,000 µL
Wash 3	D	80% Ethanol	2,000 µL
Elution plate			
Elution	A	Elution Solution	250 µL
Tip Comb	B	6-Tip Comb in 24 Deep-well plate	

3 Prepare Binding Bead Mix

- Vortex Beads vigorously to ensure they are homogenous.
- Prepare Binding Bead Mix according to the following table and sample input volume:

Component	Volume per well ^[1]
500 µL sample input	
Binding Solution	1,250 µL
Total Nucleic Acid Magnetic Beads	50 µL
Total volume	1,300 µL
500 µL – 2 mL sample input	
Binding Solution	2,700 µL
Total Nucleic Acid Magnetic Beads	100 µL
Total volume	2,800 µL

^[1] Use 10% Overage calculation when making a master mix for use with multiple samples.

- Mix well by inversion, then store at room temperature.

Digest with Proteinase K, then elute TNA

- a. Add Proteinase K to each well in Row A of the Sample Plate, according to the following table and sample input volume.

For Sample input volume	Add Proteinase K
500 µL	25 µL
>500 µL – 1 mL	50 µL
>1 mL – 2 mL	100 µL

- b. Add 500 µL to 2 mL of each sample to wells with Proteinase K in the Sample Plate.

Note: Recommend up to 1 mL input for whole blood.

- c. Invert Binding Bead Mix gently to mix, then add the appropriate amount of Binding Bead Mix to each sample in the Sample Plate according to the following table and sample input volume.

For Sample input volume	Add Binding Bead Mix
500 µL	1,300 µL
>500 µL – 2 mL	2,800 µL

Note: Remix the Binding Bead Mix by inversion frequently during pipetting to ensure even distribution of beads to all samples or wells. The mixture containing the Binding Beads is viscous. Therefore, pipet slowly to ensure that the correct amount is added. DO NOT use a repeat pipet to add to the samples as the high viscosity will cause variations in volume added.

- d. Select the program **MVP_Duo_LV** on the instrument.
- e. Start the run, then load the prepared plates into position when prompted by the instrument.
- f. After the protocol is complete (~45 minutes after start), immediately remove the elution plate from the instrument and cover the plate or transfer the eluate to a tube or plate of choice for final storage.

The purified nucleic acid is ready for immediate use. Alternatively, store the plate at –20°C for long-term storage.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Life Technologies Corporation | 2130 Woodward Street | Austin, TX 78744

For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. MAN0018073

Revision	Date	Description
B.0	06 December 2019	Updated Total Nucleic Acid Binding Buffer to Binding Solution.
A.0	14 March 2019	New document.

Important Licensing Information: This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.